


# Specification Datasheet

Viamed

SpiroTrue H

OEM Disposable Single Patient Adult Flow Sensor



CE

<b>Part numbers</b>	4310003 pack of 6 4310013 OEM pack of 20
<b>Model</b>	SpiroTrue H – OEM disposable single patient adult flow sensor + tubing set
<b>Primary packaging</b>	4310003 Box of 6 sets  4310013 PE pouch (bag) with 20 sets includes flow sensor with double PVC (DEHP free) tube connected plus adapter and hook

## Technical Data

<b>Measuring principle</b>	Differential pressure measurement
<b>Reproducibility</b>	≤ 10% at flow range 0 to 100l/min
<b>Designed in accordance with</b>	93/42/EEC (Annex I and Annex II)
<b>Housing material</b>	ABS (Acrylnitril-Butadien-Styrol)
<b>Weight</b>	≤ 9.5g (without tubing)
<b>Dead space</b>	≤ 10ml
<b>Connectors</b>	15mm O.D./15mm O.D. in accordance with ISO 5356-1
<b>Conforms to</b>	EN ISO 10993-5
<b>Packaging</b>	Non-sterile
<b>Recommended storage temperature</b>	-20°C to 50°C
<b>Humidity (R.H.)</b>	5% to 95%
<b>Storage period</b>	≤ 3 years (unopened, away from direct sunlight)
<b>Operating temperature</b>	10°C to 40°C
<b>Relative humidity</b>	5 to 95%, non-condensing

**Notes:** Before using the sensor consult the appropriate ventilator user instructions for more detailed information.

SpiroTrue H is compatible with Hamilton single use flow sensor part number 279331