



Specification Datasheet

Viamed SpiroTrue E OEM Disposable Single Patient Adult Flow Sensor



CE

Part numbers	4310065 OEM pack of 6
Model	SpiroTrue E – OEM disposable single patient adult flow sensor
Primary packaging	PE pouch (bag) with 6 pieces or 6 sets
Secondary packaging	??Box of 10 bags x 20pcs/sets each. 200pcs/200 sets total per box

Technical Data

Sensor measuring principle	Differential pressure measurement
Sensor reproducibility	≤ 10% at flow range 0 to 100l/min
Sensor housing material	ABS (Acrylnitril-Butadien-Styrol)
Sensor weight	≤ 9.5g (without tubing)
Sensor dead space	≤ 10ml
Sensor connectors	15m/15m in accordance with DIN EN ISO 5356-1
Recommended storage temperature	-20°C to 50°C
Humidity (R.H.)	5% to 95%
Recommended period of storage	≤ 3 years (unopened packaging and kept away from sunlight)
Recommended operating temperature	10°C to 40°C
Relative humidity	≤ 95%, non-condensing
Warranty period	12 months from date of sales invoice

Sensor is in non-sterile packaging

Sensor designed in accordance with 93/42/EEC (Annex I and Annex II)

Sensor conforms to EN ISO 10993-5

Sensor components including sensor tubing conforms to 2011/65/EC (RoHS)

Please note: Before using the sensor consult the appropriate ventilator user instructions for information regarding the risks related to the use of this product.

SpiroTrue E is compatible with Envitec single use flow sensor SpiroQuant H part number 07-00-0001