

Please review the data below. Once complete, click Proceed to continue.

Proceed to Next Step

Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
0022171	Universal mounting plate with `V` mount	2	✘	0.3	0.423	0.2	0.2		0	Update
0110019	Oxygen sensor R-17Vi	40	✔	0.05	0.384	0.038	0.038		1.52	Update
0110041	Oxygen sensor R-41V	61	✔	0.1	0.1125	0.039	0.039		2.379	Update
0110049	Oxygen sensor R-49V	6	✔	0.1	0.108	0.03	0.033		0.18	Update
0110057	Oxygen sensor T-7	3	✔	0.05	0.1078125	0.032	0.032		0.096	Update
0110132	Oxygen sensor - Automotive R-22AV	10	✔	0.05	0.12	0.04	0.04		0.4	Update
0110182	Oxygen sensor - Automotive R-22AVG	10	✔	0.1	0.12	0.04	0.04		0.4	Update
0110214	Oxygen sensor R-33S	200	✔	0.05	0.144	0.04	0.04		8	Update
0110215	Oxygen Sensor R-33S1	200	✔	0.05	0.1125	0.04	0.04		8	Update
0110217	Oxygen sensor R-17VAN	2390	✘	0	0	0.04	0.04		0	Update
0110220	Oxygen sensor R-22OMS	3	✘	0	0	0.04	0.04		0	Update
0110224	Oxygen sensor R-33DE	594	✘	0	0	0.04	0.04		0	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
0110225	Oxygen sensor R-22DEM	5	✘	0	0	0.04	0.04		0	Update
0110226	Oxygen sensor R-33VAN	14	✘	0	0	0.04	0.04		0	Update
0110227	Oxygen sensor R-22AEX	20	✘	0	0	0	0.07		0	Update
0110229	Oxygen sensor R-22SHO	50	✔	0.05	0.12	0.04	0.04		2	Update
0110231	Oxygen Sensor R-22DSX	10	✔	0	0	0.04	0.04		0.4	Update
0110232	Oxygen sensor R-22VAN	395	✘	0	0	0.04	0.04		0	Update
0110235	Oxygen sensor R33VAN-4	11	✔	0	0	0.04	0.04	0.080	0.44	Update
0110236	Oxygen sensor R-33S1 (RP)	10	✔	0	0	0.04	0.04		0.4	Update
0110238	Oxygen sensor R-22AT	94	✔	0	0	0.04	0.04		3.76	Update
0110239	Oxygen sensor R-33VAN (RP)	132	✔	0	0	0.04	0.04		5.28	Update
0110241	Oxygen sensor R-33VAN-2	36	✔	0	0	0.04	0.04		1.44	Update
0110243	Oxygen sensor R-22DEOX	2	✔	0	0	0.04	0.04		0.08	Update

Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
0110248	Oxygen sensor R-9507M	4	✓	0	0	0.04	0.04		0.16	Update
0110301	Oxygen sensor Type 22. P/N 3001	2550	✓	0	0	0.04	0.04	0.040	102	Update
0110381	Oxygen sensor WGA-02 Sumner	100	✓	0	0	0.04	0.04		4	Update
0110382	Sarcostopper R-17	21	✓	0	0	0.04	0.04		0.84	Update
0110383	ERFI France R-17	32	✓	0	0	0.04	0.04		1.28	Update
0110390	Oxygen sensor STK-01	100	✓	0	0	0.04	0.04		4	Update
0110420	Oxygen sensor MAX-14 (Fluke version)	7	✓	0.05	0.1125	0.038	0.038		0.266	Update
0110425	Oxygen sensor MAX-250. With 'O' ring.	17	✓	0.05	0.1	0.038	0.038		0.646	Update
0110427	Oxygen sensor MAX-250+	15	✓	0.05	0.423	0.037	0.037		0.555	Update
0110429	Oxygen sensor MAX-250E	29	✓	0.05	0.082666666666667	0.046	0.046		1.334	Update
0110560	Oxygen sensor OOM111	126	✓	0.05	0.144	0.038	0.038		4.788	Update
0110566	Oxygen sensor OOM109-LF2	1	✓	0.1	0.144	0	0		0	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
0111277	Oxygen monitor MySign O Kit	41	✓	0.8	1.5	0.33	0.33		13.53	Update
0120080	Oxygen Sensor Input Protection Cap	650	✗	0.1	0.52875	0.05	0.05		0	Update
0121350	MySign O - Oxygen Sensor Cable	11	✓	0.1	0.74025	0.5	0.5		5.5	Update
0131203	Oxygen Sensor Cable Yellow	13	✓	0	0	0.1	0.1	0.1	1.3	Update
0140009	Labour Charge - VN202 mkII	1	✓	0	0	0.01	0.01		0.01	Update
0140011	Oxygen sensor evaluation.	57	✗	0.01	1.6666666666667E-10	0	0		0	Update
0150014	Single wall cardboard boxes (305x229x229 mm) - pack of 30	10	✓	0	0	5	5		50	Update
0150030	Oxygen sensor barrier bag 3.5 x 5"	10	✗	0.1	0.144	0.01	0.01		0	Update
7910100	Oxygen Analyser - VN202 MkII Auto	200	✓	0	0	0.1	0.1	0.124	20	Update
7910110	Oxygen Analyser VN202mkII Manual	100	✓	0	0	0.1	0.1		10	Update
7910201	Oxygen Analyser Tek-Ox Auto + Quick-Ox	1000	✓	0	0	0.1	0.1	0.200	100	Update
7910210	Oxygen Analyser Tek-Ox Manual	50	✓	1	1.664	0.18	0.235		9	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
7910254	TEK-OX for Weld Gas - 0.00 UD	100	✓	0	0	0.8	0.8		80	Update

7910255	Weld Gas Analyser-Sarco Stopper 0.00 UD	200	✓	0	0	0.8	0.8		160	Update
7930007	Battery Lead	1	✓	0	0	0	0		0	Update
7930119	VN202 Mk2 LCD & Display PCB Assembly.	1	✓	0.5	0	0	0		0	Update
7930120	TEK-OX and VN202 Mk2 LCD	1050	✓	0	0	0	0		0	Update
7980000	VN202 MkII Calibration	34	✗	0.01	1.6666666666667E-10	0	0		0	Update
7990112	Vandagraph Oxygen Sensors Leaflet	100	✗	0	0	0	0		0	Update
8010100	Oxygen Sensor Sample for R&D Evaluation.	1	✓	0.1	0.12	0	0		0	Update
9070005	Protective Black Case for TEK-OX	5	✗	0	0	0	0		0	Update
9070006	Protective Case for TEK-OX	2	✓	0.1	0.144	0	0		0	Update
9071007	Jack Plug - 750 - 3.55mm	10	✗	0.1	0.144	0	0		0	Update
9071040	Molex, PicoBlade Female Connector Housing for R-9507M	50	✓	0	0	0	0		0	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
9071041	Molex PicoBlade Series Female Crimp and Wire for R-9507M	50	✓	0	0	0	0		0	Update
9530000	Label - TEK-OX Front 31x25mm	2000	✓	0.1	0.144	0	0		0	Update
9530001	Label - Tek-Ox Small Round Front	2000	✓	0.1	0.144	0	0		0	Update
9530002	TEK-OX Base Seal	2486	✗	0	0	0	0		0	Update
9530003	TEK-OX Main Body insert black (internal chassis)	2673	✗	0	0	0	0		0	Update
9530005	TEK-OX Black Bottom Cover	2504	✗	0	0	0	0		0	Update
9530006	TEK-OX Front Gasket	2649	✗	0	0	0	0		0	Update
9530019	TEK-OX Box - model BP-3003 - Top Label	1000	✓	0	0	0	0		0	Update
9530020	Label Cutter - One off Cost	1	✓	0	0	0	0		0	Update
9530026	TEK-OX PCB V0.6 with the Display fitted	4	✓	0	0	0	0		0	Update
9530027	VN202 PCB V0.8 with the display fitted	1	✓	0	0	0	0		0	Update
9530052	TEK-OX Left Shell Black (strengthened)	3091	✗	0	0	0	0		0	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
9530053	TEK-OX Right Shell Black (strengthened)	2502	✗	0	0	0	0		0	Update
9530055	TEK-OX Front Cover Black with Lens and Protect	2500	✗	0	0	0	0		0	Update
9711000	Flow Diverter	115	✗	0	0	0	0		0	Update

9711770	Oxygen Analyser VN202 - Service Exchange	2	✓	0	0	0	0		0	Update
9711901	DIN Pressure Reducer	40	✗	1	0.423	0	0		0	Update
9711905	Tee Blue 2 x 1 Way Valve	13	✗	0	0	0	0		0	Update
9711912	Vandagraph TEE 2 x One Way Valve	5000	✓	0	0	0	0		0	Update
9730001	Acrylic Conformal Coating 200ml	1	✓	0	0	0	0		0	Update
9730036	Push Switch Round Waterproof	1000	✓	0	0	0.01	0		10	Update
9730040	10K Calibration Pot GPB 16NP10KMA	1000	✓	0	0	0.01	0		10	Update
9730063	TEK-OX Battery Bottom Contact, Nickel plated	1823	✓	0	0	0	0		0	Update
9730064	TEK-OX Battery Terminals, left Nickel plated	1963	✓	0	0	0.01	0		19.63	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
9730067	Tek-Ox Battery Terminal Coils Right	1767	✓	0	0	0	0		0	Update
9730069	Label - TEK-OX Box 77 x 41mm	1000	✗	0.5	0.6345	0	0		0	Update
9730087	VN202 Protective Carry Case	511	✓	0	0	0	0		0	Update
9730088	Seal Kit 0098-92015	1000	✗	0.1	0.144	0	0		0	Update
9730181	R-22AT Adapter Cap	322	✓	0	0	0	0		0	Update
9730210	Quick-Ox Sampler	110	✗	0.1	0.144	0	0		0	Update
9730215	Gas Sampling DIN KIT	87	✗	1	0.6345	0	0		0	Update
9790000	Book - Practice of O2 Measurement for	1	✗	0.5	0.423	0	0		0	Update
9790002	Book - Oxygen Measurement for Divers	1	✗	0	0	0	0		0	Update
9790003	Amazon Digital Book sales	17	✗	0	0	0	0		0	Update
9790065	Tek-Ox Screw M3 x 10mm	10	✗	0	0	0	0		0	Update
9910200	Sample supplied for evaluation purposes.	1	✗	0.01	1.6666666666667E-10	0.01	0.01		0	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
9950051	Batteries 1.5V AA. Pk10.	340	✗	1.5	0.52875	0	0		0	Update
EXW	Delivery: EXW - VST, UK (Incoterms 2020)	61	✗	0	0	0.01	0.01		0	Update
PPRD	Royal Mail - 1st Class Signed For	256	✗	0.01	1.6666666666667E-10	0	0		0	Update
PPRM1C	Royal Mail 1st Class	3	✗	0.01	1.6666666666667E-10	0	0		0	Update
PPRMAS	Royal Mail Airmail - Standard	1	✗	0	0	0	0		0	Update
PPRMITS	Royal Mail International Tracked and Signed	142	✗	0.01	1.6666666666667E-10	0	0		0	Update

PPSD	Royal Mail - Special Delivery	4	✘	0.01	1.6666666666667E-10	0	0		0	Update
PPUPS1	UPS Courier Delivery - Standard	167	✘	0	0	0.01	0.01		0	Update
PPUPS2	UPS Courier Delivery - Standard	1	✘	0.01	1.6666666666667E-10	0	0		0	Update
PPUPS4	UPS Courier Delivery - Express	1	✘	0.01	1.6666666666667E-10	0	0		0	Update
PPUPS6	UPS Courier delivery - Standard.	291	✘	0	0	0.01	0.01		0	Update
PPUPS7	UPS Courier Delivery - Express Saver	69	✘	0	0	0.01	0.01		0	Update
Stock Ref	Description	Item Count	WEE Applicable	Package Weight (kg)	Volumetric Weight (kg)	WEE Weight (kg)	Physical Item Weight (kg)	Pen Wee Weight	Feedback	Action
PPUPS8	UPS Courier delivery - Express	1	✘	0	0	0.01	0.01		0	Update
REPCUSTOMER	O2 Analyser Repair	2	✘	0	0	0	0		0	Update
REPHEO2	Helium/O2 Analyser Repair	2	✘	0	0	0	0		0	Update
REPTEK	Repair TEK-OX	22	✘	0	0	0	0		0	Update
REPVAN	Repair VN202	15	✘	0	0	0	0		0	Update
Total Feedback:									643.614	